



BEEKLEY LOCALIZERS


Specific Localization | Specific Communication | Specific for CT

The first skin markers designed specifically for use in diagnostic CT! Save time and reduce radiation dose to patients by eliminating repeat exams because of skin marker artifact!

Localizers provide CT-Scan technologists with a clear and consistent means of communicating areas of interest to the interpreting radiologist.

Localizers are made of a unique and proprietary formulation that was specifically calibrated for use in diagnostic CT. This formulation enables a Localizer to image brightly on the CT scan but without potentially obscuring artifact. Localizers can increase patient throughput by helping to reduce the need to re-scan while also serving to limit radiation dose for the patient.

Localizers Image Brightly without Obscuring Artifact



Localizer specifies an exact area of interest without troublesome artifact

Older methods generate considerable artifact



Artifact from older methods can potentially lead to repeat scans or obscured areas of concern



Localizer

order code: 210

contents per box: 40

- 2.3mm low density pellet
- specific to diagnostic CT
- perfect for identifying soft tissue masses on the neck



Localizer

order code: 211

contents per box: 40

- 4.0mm low density pellet
- specific to diagnostic CT
- perfect for identifying areas of interest on abdominal cases or with larger patients

Pricing:

To begin using Localizers with your Diagnostic CT exams contact your Beekley Sales Specialist today at 1-800-BEEKLEY (233-5539) or on the web at www.beekley.com

12+ boxes
\$39.00

6-11 boxes
\$44.00

1-5 boxes
\$53.00