



Marking Specific Points and Lines of Reference on Treatment Sites

Conventional Simulation



BEEKLEY SPOTS® Metallic Markers: designed for use in conventional simulation and on port films. Lead-free metallic radiopaque pellet and linear markers are ideal for a variety of treatment planning points of reference.



T-SPOTS .8mm wire used to identify size of abdominal tumor.

- Lead-free, metallic pellet and linear markers image brightly
- A variety of different sizes in both pellet and linear shapes
- Identify and differentiate points of reference on treatment sites
 - > X-SPOTS 1.5mm, Y-SPOTS 2.0mm, V-SPOTS 2.5mm pellets
 - > Z-SPOTS 4.3mm burn-out resistant pellet
 - > T-SPOTS .8mm continuous line
 - > S-SPOTS .33mm continuous line
 - > D-SPOTS radiopaque arrow

Perfect for Marking:

- ✓ Isocenters
- ✓ 3 point set-ups
- ✓ Outer canthus, nodes
- ✓ Scars
- ✓ Field borders, breast tangents
- ✓ Match lines

X-SPOTS® Metallic Pellet




1.5mm pellets used to mark 3 point set-up on breast patient.

Y-SPOTS® Metallic Pellet



2.0mm pellets used to mark pelvic anterior and laterals for 3 point set-up on prostate patient.

T-SPOTS® and S-SPOTS® Metallic Lines



.8mm and .33mm radiopaque continuous lines used to differentiate medial and lateral breast tangents.

BEEKLEY SPOTS® Metallic Skin Markers

Designed for use in conventional simulation and on port films. Lead-free metallic radiopaque pellet and linear markers are ideal for a variety of treatment planning points of reference.



X-SPOTS® Pellet
order code: **101**
size: **1.5mm diameter**
contents per box: **150**



Z-SPOTS® Pellet
order code: **114**
size: **4.3mm diameter**
contents per box: **40**




T-SPOTS® Line
order code: **111**
size: **.8mm diameter**
contents per box: **275cm**



Y-SPOTS® Pellet
order code: **102**
size: **2.0mm diameter**
contents per box: **132**



D-SPOTS® Arrow
order code: **104**
contents per box: **99**



S-SPOTS® Line
order code: **107**
size: **.33mm diameter**
contents per box: **363cm**



V-SPOTS® Pellet
order code: **103**
size: **2.5mm diameter**
contents per box: **130**

Product shown actual size.

**SPOTS Metallic
Skin Markers**

1-4 boxes
\$98.00

5-9 boxes
\$84.00

10-19 boxes
\$80.00

20-49 boxes
\$70.00

50-99 boxes
\$53.50

100+ boxes
\$37.00

Combine & Save! Combine quantities on all types of Beekley SPOTS products to earn maximum discounts.

**For more information or to place an order
contact your Beekley Radiation Oncology Sales Specialist today at
1-800-BEEKLEY (233-5539)**