



NEW!

BEEKLEY *Tender Touch*TM Thermoplastic Edge Tape



Silky, soft, smooth, cushioned tape protects patients from the rough edges on thermoplastic masks. TenderTouch, developed specifically for this application, is perfect for those times when masks are cut for comfort, claustrophobia, bolus, bite blocks and electron fields.

New Method TenderTouchTM



Maximum Patient Comfort

Silky, soft 1" wide material covers and cushions the sharp, ragged edges of thermoplastic masks. The tape is cushy and pliable providing ultimate patient protection and comfort.

Quick & Easy to Apply

Scored edges mold easily around corners and allow for larger piece applications, saving valuable time. Dispenses conveniently from plastic package.

Stays on for Duration of Treatment

Sticks consistently to the thermoplastic material for the entire treatment period from 4 to 8 weeks.

Durable Marking Surface

For marking crosshairs, field borders, leveling lines, shielding marks, etc. that will not fall off during treatment.

Dosimetric Properties of TenderTouchTM

The effects on surface dose of the TenderTouch foam adhesive material is minimal and could be ignored in most clinical situations. In particular, the incremental change in surface dose for this material relative to the effect of the thermoplastic is not considered clinically significant in most foreseeable situations. A copy of the full Dosimetric Report is available upon request.

Old Methods Paper, Silk, & Masking Tapes



Rough Ragged Edges Protrude

Rough edges may poke through and cause patient discomfort because the paper, masking, scotch and cloth tapes are too thin and not cushioned.

Time Consuming Application

It takes valuable time to cut small sections for the curved cut out areas on the masks. Also need to apply several layers to achieve an adequate level of cushioning for patient comfort.

Tape Falls Off

Tape does not stick as well. It can fall off as it loses its adhesive over time during the treatment period and may need to be re-applied. This can be a critical accuracy issue if used to mark treatment planning reference points.

We have a soft spot for your patients!

Order Code: 418
Contents: 310cm per box
Size: 2.5cm/1" wide

Pricing:

1-4 boxes \$60/Box | 5-9 boxes \$55/Box | 10+ boxes \$50/Box

**FOR MORE INFORMATION OR TO PLACE AN ORDER,
CONTACT YOUR BEEKLEY RADIATION ONCOLOGY SALES SPECIALIST AT**

1-800-233-5539



"patent pending"