

# PINPOINT®

## Multi-Modality Fiducial Markers for CT / MRI Image Registration

fiducial markers with a unique conical design  
to assist in accurate image registration

### PinPoint multi-modality fiducial marker provides...

#### ACCURACY

Unique conical design direct probes to 1.27mm diameter center hole—PinPoint's registration area is 37% smaller than other commonly used fiducial markers.

#### UNMISTAKABLE IMAGING

Bright and consistent imaging in both CT and MRI—PinPoint is easily identified on scans.

#### TIME SAVINGS

Easy application and removal—PinPoint's Pik-Up™ tab helps eliminate frustration and wasted time.



REF 128 QTY 30/Box

\* product shown not actual size and not to scale.

**Call 1.800.233.5539 • Fax 1.800.735.1234**  
**Visit [beekley.com](http://beekley.com) • Email [info@beekley.com](mailto:info@beekley.com)**

Customers outside the U.S. – Contact your local distributor for pricing and product availability.  
To locate a distributor call +1.860.583.4700 or email [international@beekley.com](mailto:international@beekley.com)

Manufactured by Beekley Corporation  
One Prestige Lane, Bristol, CT 06010-7454 USA  
Tel: 1.800.233.5539 or +1.860.583.4700  
Fax: 1.800.735.1234 [beekley.com](http://beekley.com)

BEEKLEY, BEEKLEY MEDICAL, PINPOINT, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada. PIK-UP is an unregistered trademark of Beekley Corporation.

© 2023 Beekley Corporation. All rights reserved.

#### WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Not made with natural rubber latex
- Lead-Free
- Non-Sterile



EMERGO EUROPE  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands



REV: PP128\_INF\_0723