



for use with IONIC contrast
(contrast to be added separately)
Tropical Fruit Flavor, Sugar-Free



Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010 USA
Tel: 1-800-233-5539 or +1-860-583-4700 Fax: 1-800-735-1234 www.beekley.com

BEEKLEY, BEEKLEY MEDICAL, BREEZA, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada. All other trademarks are the property of their respective owners. Pat. - www.beekley.com/patents © 2014 Beekley Corporation. All rights reserved.

Description

Breeza for use with IONIC oral iodinated contrast is a flavored beverage used in CT-Scan.

REF 220

16.9 FL OZ (473 mL)

Indications

Breeza for use with IONIC oral iodinated contrast is a flavored beverage used in CT-Scan.

Warnings and Precautions

- Allergens: may contain traces of fruit
- Do not drink if allergic to any of the ingredients
- Use only under the supervision and direction of a physician
- Discard any unused portion after opening
- Do not use if tamper-proof seal is separated prior to opening

User Guidance

- For best results, refrigerate prior to using
- Use immediately after adding contrast
- Use by expiration date stamped on bottle cap

Instructions for Use

1. Fill out patient information section either directly on Breeza bottle or on enclosed labels using a fine tip Sharpie® marker. Adhere label to bottle.
2. Twist cap open and add iodinated contrast as per recommended protocol.
3. Firmly twist cap back on; shake gently to mix.
4. Repeat steps 1–3 for any additional bottles as necessary.
5. Deliver prepared bottle(s) of Breeza to patient.

Ingredient Information

NUTRITION FACTS

Serving Size: 8 FL OZ (237 mL)
Servings per Container: 2

AMOUNT PER SERVING

Calories: 0

Calories from Fat: 0

	% Daily Value*
Total Fat 0g	0%
Sodium 120mg	5%
Potassium 10mg	<1%
Cholesterol 0mg	0%
Total Carbohydrates <1g	<1%
Sugars 0g	
Sugar Alcohol <1g	
Protein 0g	0%

(Not a significant source of calories from fat, saturated fat, trans fat, cholesterol, dietary fibers, vitamin A, vitamin C, calcium, or iron.)

*Percent Daily Values are based on a 2,000 calorie diet.

Ingredients: Filtered Water, Sorbitol, Citric Acid, Natural and Artificial Flavors, Sodium Hexametaphosphate (preserves freshness), Mannitol, Sodium Citrate, Malic Acid, Sodium Benzoate and Potassium Sorbate (preserves freshness), Sucralose, Acesulfame K, Calcium Disodium EDTA.