



The Beekley Skin Marking System[®] for Mammography


A visual communication system of radiographically distinct skin markers.




Pellet = Nipple




Circle = Mole



Line = Surgical Scar



Triangle = Palpable Mass



Square = Area of Concern or Pain



beekley.com • 800.233.5539

Visit beekley.com for Product Safety Information

Beekley Medical[®]
when your diagnosis must be right[®]



Manufactured by Beekley Corporation, Bristol, CT 06010-7454 USA

BEEKLEY, BEEKLEY MEDICAL, BEEKLEY SKIN MARKING SYSTEM, WHEN YOUR DIAGNOSIS MUST BE RIGHT, the Raised Circle, the Raised Triangle, the Raised Square, and the Hummingbird design are Reg. U.S. Pat. & Tm. Off. BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada.
© 2023 Beekley Corporation. All rights reserved.