


Ortho-SPOT[®]


BEEKLEY MEDICAL[®]
when your diagnosis must be right[®]

MRI Skin Markers for Orthopedic Imaging

Ortho-SPOT Packets[™] minimize indentation and the risk of distorting pathology; clearly communicate area of interest:



Capsule markers on extremities create indentation and can potentially obscure pathology, raising questions regarding exact area of interest:



.75cm Radiance[®] filled Ortho-SPOT Packets[™]

REF 186 QTY 40/Box

PACKET MARKERS ONLY
FROM BEEKLEY MEDICAL

- Precise, indentation-free localization
- Sized to image only on slices specific to area of interest
- Bright consistent imaging on every sequence



Call 1-800-233-5539 • Fax 1-800-735-1234
Visit www.beekley.com • Email info@beekley.com

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1-860-583-4700 or email international@beekley.com

Manufactured by Beekley Corporation
One Prestige Lane, Bristol, CT 06010 USA
Tel: 1-800-233-5539 or +1-860-583-4700
Fax: 1-800-735-1234 www.beekley.com

BEEKLEY, BEEKLEY MEDICAL, SPOT, ORTHO-SPOT, RADIANCE, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off. BEEKLEY and SPOT are Registered Community Trademarks. BEEKLEY is a registered trademark in Canada. ORTHO-SPOT PACKETS, the PACKET product configuration, and PIK-UP are unregistered trademarks of Beekley Corporation.
© 2017 Beekley Corporation. All rights reserved.

WARNINGS AND PRECAUTIONS:

- Single Patient Use Only; Do Not Reuse 
- External Use Only; Do Not Ingest or Swallow
- Do Not Place on Compromised Skin
- Remove After Use
- Latex-Free
- Lead-Free
- Non-Sterile 

CE EC REP EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

