

PINPOINT®

Multi-Modality Fiducial Markers for CT / MRI Image Registration

breakthrough fiducial markers with a unique conical design to assist in accurate image registration

PinPoint provides...

ACCURACY

Unique conical design direct probes to 1.27mm diameter center hole—PinPoint's registration area is 37% smaller than other commonly used fiducial markers.

UNMISTAKABLE IMAGING

Bright and consistent imaging in both CT and MRI—PinPoint is easily identified on scans.

TIME SAVINGS

Easy application and removal—PinPoint's Pik-Up tab helps eliminate frustration and wasted time.



REF 128 QTY 30/Box

* product shown not actual size and not to scale.

Call 1-800-233-5539 • Fax 1-800-735-1234
Visit www.beekley.com • Email info@beekley.com

Customers outside the U.S. – Contact your local distributor for pricing and product availability.
To locate a distributor call +1-860-583-4700 or email international@beekley.com

BEEKLEY MEDICAL®
when your diagnosis must be right®

One Prestige Lane, Bristol, CT 06010 Tel: 1-800-233-5539 or +1-860-583-4700 Fax: 1-800-735-1234 www.beekley.com

BEEKLEY, BEEKLEY MEDICAL, PINPOINT, and WHEN YOUR DIAGNOSIS MUST BE RIGHT are Reg. U.S. Pat. & Tm. Off.
BEEKLEY is a Registered Community Trademark. BEEKLEY is a registered trademark in Canada.

© 2014 Beekley Corporation. All rights reserved.



REV: PP128_INF_1014